

## A. NeuroChaT GUI Use Guide

### Download and installation:


Download the codes from OSF website or GitHub

- OSF link: <https://osf.io/kqz8b>
- GitHub link: <https://github.com/mnislamraju/NeuroChaT>

### Start NeuroChaT GUI (no installation required)

- Using command prompt:
  - Change your current directory to the NeuroChaT directory
    - >cd C:\path\to\neurochat
  - Execute following command
    - >python neurochat.py
  - Or, >python C:\path\to\neurochat\neurochat.py
- Using Python IDE, i.e. Spyder:
  - Open neurochat.py in the editor
  - Press the **Run** button to start NeuroChaT GUI

Following window will pop-up:




*Figure 1: NeuroChaT graphical user interface*

## Overview of the front panel:

### Input Data Format

Select one of the three data or file formats. NeuroChaT is currently supporting Axona, Neuralynx and NWB (HDF5) file formats.



*Figure 2: Select the input file format from one of the three options in the 'Input Data Format' dropdown menu.*

## Analysis Mode

NeuroChaT works in three analysis modes (**Figure 3**): Single Unit, Single Session, Listed Units. Following are the descriptions of how these modes work.

**a. Single Unit:** When selected, it analyses data that belongs to a single cluster extracted from spike-sorting methods. This mode is a perfect choice when the detailed analysis using many functions is required to explore more properties of the cell.

**b. Single Session:** This mode analyses data from a single spike file. In Neuralynx system, it accepts .ntt and .nst files. For Axona system, it supports .n (n= 1, 2, 3 etc.) files. This mode looks for all the clusters that have been identified from the spike sorting methods, and analyse the selected functions for all of them. Along with appropriate choice of analysis functions, this is a powerful mode for thorough examination and characterization of units from a single electrode.

**c. Listed Units:** This mode analyses data from units those are listed in an Excel format as shown in **Figure 3**. The first column is the directory where the data is stored. Data specifications are provided by following means:

**HDF5 File:**

Column 1	Column 2	Column 3	Column 4
Data directory	Name of the HDF5 file without extension	Single unit of interest	LFP channel ID

#### Axona and Neuralynx Files:

Column 1	Column 2	Column 3	Column 4	Column 5
Data directory	Name of the spatial data file without extension	Name of the spike data file with extension	Single unit of interest	LFP channel extension (Axona)  Or Name of the LFP data file (Neuralynx)

	A	B	C	D
1	dir	nwb	cell id	lfp
2	C:\Users\ [REDACTED]	120213_26	6	eeg


	A	B	C	D	E
1	Directory	position file	Spike file	cell id	lfp chan
2	C:\Users\ [REDACTED]	120213_26_1	120213_26.6	4	eeg

**Figure 3:** Input style in Excel files for batch mode analysis using 'Listed Cell'. Top row shows the style for HDF5 data. Bottom row shows the style for Axona and Neuralynx systems.

Note: For Axona system, you should provide the spatial information as a .txt file format. If you use TINT for spike-sorting, you can export these information from there as .txt file.


#### Unit No

Select the unit that you are interested to analyse. Although NeuroChaT lists up to 255 units, units those are identified using the clustering process can only be used.



## LFP Ch No

Select the LFP channel from the dropdown list. When the data is browsed using the 'Browse' button this box fills with the potential LFP channels in the folder with recorded data. For HDF5 dataset, it shows the data groups in the directory '/processing/Neural Continuous/LFP'. It is always 'ncs' For Neuralynx data to refer to the files with .ncs extension.



*Figure 4: Dropdown list to select the LFP channel*

## Graphic Format

Select one of the two options for the file format of graphical output from NeuroChaT analyses.

## Browse

Clicking on this button prompts for the file or folder selection based on analysis mode. For 'Single Unit' and 'Single Session' mode, this will ask user to select .ntt/.nst file followed by .nvt file (Neuralynx) or .n file followed by .txt file (Axona). For 'Listed Units' mode, it will ask for specifying the .xls/.xlsx file that contains the list of the units.

## Analysis Selection

This section provides a list of analyses that can be selected by ticking the boxes beside their names. Checking the 'Select All' box selects all the analyses and unchecking it removes their selection.

## Select Cell Type

Analyses of interest can also be set by pressing one of the buttons in the ‘Select Cell Type’ section. It selects the analyses of interests essential for individual unit type. For example, characterizing units for rhythmic properties does not require spatial analyses. Therefore, only the analysis set that characterises spike trains are selected. The user can select or deselect analyses from this set.

### **Log Box**

The box with white background at the bottom of the user-interface displays log of NeuroChaT actions, warnings and errors. Warnings are represented by orange texts while errors are displayed in red texts. All other NeuroChaT logs are shown as blue texts.

### **Save Log**

Pressing ‘Save Log’ button prompts the user to save the texts at the log box. This will export the log texts as plain ASCII texts and, therefore, there will be no colour in the output text file.

### **Clear Log**

This button simply clears the log record in the log box. It is recommended that the log box is cleared at intervals so that errors are easier to find when the logs are exported

### **Start**

Pressing this button starts reading the data from the specified files, and running analyses on them. As long as the data analysis keeps going, this button remains disabled to avoid unnecessary interruption of the execution. It is enabled again after the completion of the analysis or if there is a fatal error which stops the execution. At the end of the analysis of each unit, log box shows the full file directory of the output graphics which are saved in native folder of the spike data. If the selected functions have numeric results, a table appears after successful execution of all analysis functions. A sample table is shown in **Figure 5**. Clicking on the ‘Export results’ prompts the user to save the tabular data in Excel file (**Figure 5, Bottom**).

	Mean Spiking Freq	Mean width	Mean amplitude	Mean height	Std height	Std width
TT6_SS_4_eeg	9.73022481266	241.153944891	203.199721643	204.495651026	21.4653091143	64.52079776
TT3_SS_1_eeg	21.0133222315	126.49912166	241.516174434	241.710288474	18.5259256623	13.3141324968
TT7_SS_2_eeg	0.706666666667	358.195754717	172.367887677	182.170939895	27.0331326065	48.0785555413
TT5_SS_1_eeg	35.682764363	189.698790495	357.959282513	743.145034856	25.8458682168	15.9444303437

	A	B	C	D	E
1		Mean Spiking Freq	Mean width	Mean amplitude	Mean height
2	TT6_SS_4_eeg	9.730224813	241.1539449	203.1997216	204.495651
3	TT3_SS_1_eeg	21.01332223	126.4991217	241.5161744	241.7102885
4	TT7_SS_2_eeg	0.706666667	358.1957547	172.3678877	182.1709399
5	TT5_SS_1_eeg	35.68276436	189.6987905	357.9592825	743.1450349
6	TT4_SS_1_eeg	54.81681932	144.4640642	148.4015041	323.181652

Figure 5: Sample of parametric results of NeuroChaT analyses in a tabular form. The table pops up once the analyses are complete. Bottom row shows how it looks like once it is exported to an Excel file.

## Menu items

Current version of NC has 4 main menus items: File, Settings, Results & Help.

### File

File menu consists of *Open*, *Save Session*, *Load Session* and *Exit* options as shown in **Figure 6** below.

**a. Open:** This acts exactly as browse button described before. It will prompt user to give appropriate input depending on the analysis mode selected. If no file is selected, it generates a warning.

**b. Save Session:** This will allow the user to save the configuration of NeuroChaT in *.ncfg* file (i.e. input and output format, analysis mode, cell no, file and directory information from user input, the state of the functions that are selected for analysis, and the parameters used for the analyses). If aborted, it creates a warning in the log box.

**c. Load Session:** This option will take the user to load configuration (*.ncfg*) file. If no file is selected, it will show a warning in the log box.

**d. Exit:** Terminates the NeuroChaT software.






Figure 6: 'File' menu items. Shortcuts for invoking each of these actions are also shown.

## Settings

Settings menu has one item: *Parameters*. Clicking this item initiates the parameter selection box as shown in **Figure 7**. Analysis specific parameters and their accepted values are listed in **Appendix D**. Clicking on an item on the left panel will display the parameter setting panel for that analysis on the right panel.




Figure 7: Window for analysis specific input parameter settings.

## Utilities

This menu comprise of useful NC utilities described below (Figure 8):




Figure 8: Items in the 'Utilities' menu.


**a. Export results:** This action prompts the user to save recent analysis results in Excel file.


**b. Export graphics file info:** The user can export the directory and name of the output PDF/PS files containing analysis graphics and the HDF5 files associated to the data in an Excel file. The sample output looks as follows:

	A	B	C
1		<b>Graphics Files</b>	<b>NWB Files</b>
2	<a href="#">TT6_SS_4_eeg</a>	<a href="#">\120213_26_TT6_SS_4_eeg.pdf</a>	<a href="#">\120213_26.hdf5</a>
3	<a href="#">TT3_SS_1_eeg</a>	<a href="#">\120412_1_TT3_SS_1_eeg.pdf</a>	<a href="#">\120412_1.hdf5</a>
4	<a href="#">TT7_SS_2_eeg</a>	<a href="#">\062612_1_TT7_SS_2_eeg.pdf</a>	<a href="#">\062612_1.hdf5</a>
5	<a href="#">TT5_SS_1_eeg</a>	<a href="#">\112512_1_TT5_SS_1_eeg.pdf</a>	<a href="#">\112512_1.hdf5</a>
6	<a href="#">TT4_SS_1_eeg</a>	<a href="#">\102612_1_TT4_SS_1_eeg.pdf</a>	<a href="#">\102612_1.hdf5</a>

*Figure 9: Sample Excel file from the export of graphics file information.*

**c. Merge output PDF/PS:** This utility allows merging multiple PDF/PS files into a single PDF or PS file. The utility works in two ways, as shown on the top row of **Figure 10**: 1. By uploading an Excel list of the PDF/PS file names (full name, including directory) using the ‘Use Excel List’ option, 2. By picking the files manually using ‘Choose files’ option. ‘Select now’ button will be activated, and clicking on this will display a file-picking utility. The ‘Save In’ button opens a save file dialogue for the user to select the file. The utility executes upon pressing the ‘Start’ at the bottom of its window. Origin files remain intact.





*Figure 10: Upper row: Merge file utility showing the options of 'Use Excel List' or 'Choose Files'. Bottom Row: The file-picker for selecting files to merge*

**d. Accumulate output PDF/PS:** This utility help accumulating the analysed output graphics files into a single folder. This works in the same way as merging files described above i.e., takes Excel file input or allows manual pick up. Specified files are then accumulated in the folder selected by 'Save In' button (**Figure 10, bottom**). This utility executes when the 'Start' button, located at the bottom of the window, is pressed. PDF/PS files being accumulated are not deleted from their original location.

**e. Verify units:** Clicking this item asks the user to upload an Excel file that contains the name of the directory, spike file, and the unit number that the user is verifying as shown by the sample in **Figure II**. Output of this analysis adds the last two columns that shows if the file and the unit in that file exists.

	A	B	C	D	E	F
1		Directory	Spike file	cell id	fileExists	unitExists
2	0		120213_26.6	4	TRUE	TRUE
3	1		120412_1.3	1	TRUE	TRUE
4	2		062612_1.7	2	TRUE	TRUE
5	3		112512_1.5	1	TRUE	TRUE
6	4		102612_1.4	1	TRUE	TRUE

*Figure II: Sample input format and output in 'Verify units' utility.*

**f. Evaluate clustering:** This option evaluates the quality of clustering by measuring Bhattacharyya coefficient and Hellinger distance between the clusters of a recording session. The clusters are formed using peaks, troughs and two principle components in each electrode channels. Sample input and results are shown in **Figure 12**.

	A	B	C	D	E	F
1		Directory	Spike file	cell id	BC	Dh
2	0		120412_1.3	1	1.16E-12	1
3	1		062612_1.7	2	5.04E-08	1
4	2		112512_1.5	1	8.72E-19	1

**Figure 12:** Sample Excel file input for cluster evaluation. The analysis output are added at the tail of each row of data.

**g. Compare single units:** This option allows the user to compare the units in two different recordings. The analysis takes the clusters, formed as that of cluster evaluation, from two sessions, and compares their similarity by measuring their degree of overlap (Bhattacharyya coefficient and  $\chi^2$ -similarity) or their statistical distance (Hellinger distance). Sample input style and results are shown in **Figure 13**.

	A	B	C	D	E	F	G	H	I
1	Directory 1	Spike File 1	Cell 1	Directory 2	Spike File 2	Cell 2	BhattCoeff	Hellinger	Chi Sq Sim
2		061516_ELSSW1.4	1		061716_ELSNW2.4	1	0.393037264	0.779078	0.31384119
3		061516_ELSSW1.4	1		061816_ELSW3.4	1	0.430171872	0.75487	0.32798797
4		061516_ELSSW1.4	1		061516_ELSSW1.4	1	1	0	1

**Figure 13:** Sample Excel file input for cluster similarity measurement. Two sets of specifications are required for comparison. The analysis output are added at the tails of each set of specifications of the comparing units.

**h. Convert to NWB formats:** This option takes a list of file information in Excel format, as shown in **Figure 14**, and converts the data into HDF5 (NWB) format. The sample is for Axona recordings. Full filename (without directory) should be written in LFP data specification columns.

	A	B	C	D
1	Directory	Behavioural file	Spike file	lfp chan or file
2		120213_26_1	120213_26.6	eeg
3		120412_1_1	120412_1.3	eeg
4		062612_1_1	062612_1.7	eeg

**Figure 14** Sample Excel input for converting the recorded data to HDF5 format.

In all of these Excel based manipulations, columns names are not strictly defined, but the order of information is strictly followed. The first row is always considered to specify the header for the Excel data.

## B. API Use Guide

This API use guide produces results reported in **Chapter 5**.

In addition to the codes for verifying the place cell, head-directional cell and analyses of rhythmic units, it also shows examples of other useful methods that can be harnessed for creating simple and efficient analysis scripts and data management.

Please refer to the code-documentation for the description of each module, their classes and functions, and methods in each class.

In addition to the example units, this guide shows uses of NeuroChaT and its components in many different ways.

### Step-1: Download NeuroChaT package from [OSF](#) or [GitHub](#)

NeuroChaT can be used without any burden of installation. You can download a local copy of NeuroChaT codes and insert the path to your system.

### Step-2 Insert NeuroChaT path to \$PYTHONPATH

```
import sys
sys.path.insert(1, 'path\to\neurochat')
```

### Step-3 Import modules and classes

We are importing only NSpike and NSpatial for the moment. We will add and import NLfp data for analyses that require LFP signals. nc\_plot is the module that provides with plotting functions

```
from neurochat.nc_data import NData
from neurochat.nc_spike import NSpike
from neurochat.nc_spatial import NSpatial
import neurochat.nc_plot as nc_plot
```

## Step-4 Instantiate objects

The names C0 and S0 for the unit and the spatial data are arbitrary

```
spike= NSpike(system = 'Axona')
spike.set_name('C0')

spat= NSpatial(system = 'Axona')
spat.set_name('S0')
```

## Step-5 Add names for the data files

```
data_dir= '\full\file\directory\of\place cell\recorded\Axona\data\'

spat.set_filename(data_dir + '040513_1_1.txt')
spike.set_filename(data_dir + '040513_1.6')
```

### For HDF5 files,

Path of the data should also be added following a '+' sign. The system argument should be changed or could be set at NSpatial(system= 'Axona')

```
spat.set_system('NWB')
spike.set_system('NWB')

data_dir= '\full\file\directory\of\place\cell\HDF5\data\'
spat.set_filename(data_dir + '040513_1.hdf5+/processing/Behavioural/Position')
spike.set_filename(data_dir + '040513_1.hdf5+/processing/Shank/6')
```

## Step-6 Load spatial and spike data. Set the unit number

```
spat.load()
spike.load()

spike.set_unit_no(3)
```

## Step-7 Instantiate NData object. Add individual data objects to NData object.

```
ndata= NData()
ndata.spike= spike
ndata.spatial= spat
```



The data format, filenames for individual datasets can be set using `ndata`

```
ndata.set_data_format(data_format = 'NWB')  
  
ndata.set_spatial_file(data_dir + '040513_1.hdf5+/processing/Behavioural  
/Position')  
ndata.set_spike_file(data_dir + '040513_1.hdf5+/processing/Shank/6')
```

They can be loaded using `ndata`

```
ndata.load()
```

Or, individually

```
ndata.load_spatial()  
ndata.load_spike()
```

And the unit number can be set as well

```
ndata.set_unit_no(3)
```

## Step-8 Perform analysis of interest

### Analysis of place cell

**Place cell firing map by using `ndata`:**

Pixel size is set 3cm. A 5x5 box filter is used for smoothing the firing map

```
placeData= ndata.place(pixel = 3, filter = ['b', 5])
```

**Similar results can be obtained by passing timestamps of the spiking unit to the `spat.place()` method**

`NData` object performs the job of connecting these two objects and simplifies the analysis


```
placeData= spat.place(spike.get_unit_stamp(), pixel = 3, filter = ['b',  
5])
```

### Plotting relevant data

Refer to `neurochat.nc_plot.py` module for more plotting functions.

*Following command is used for inline display of graphics in Notebook*

```
%matplotlib inline  
fig= nc_plot.loc_firing(placeData)
```




### **Analysis and plotting of locational shuffling analysis using default parameters**

The spike timestamps are shuffled for 500 times. Pixel size is 3 cm. limit=0 implies that the spikes timestamps are randomly shuffled in (-duration, +duration) range

```
pshuffleData= ndata.loc_shuffle(nshuff = 500, limit = 0, pixel = 3)  
fig= nc_plot.loc_shuffle(pshuffleData)
```


## Distribution of locational firing specificity indices




**Analysis and plotting of locational shifting analysis using shifting index from -10 to +20**

Spike timestamps are gradually shifted from -10 to +20 units of spatial time-resolution. If the video for tracking animal behaviour is sampled at 50Hz, this means the spike-train is shifted from -200ms to +400ms


```
import numpy as np # numpy imported for the use of np.range
pshiftData= ndata.loc_shift(shift_ind = np.arange(-10, 20))
fig= nc_plot.loc_time_shift(pshiftData)
```






### Head directional analysis of this unit

```
hdData= ndata.hd_rate()
fig= nc_plot.hd_firing(hdData)
```



### Multiple regression analysis and plotting

```
regressData= ndata.multiple_regression()
fig= nc_plot.multiple_regression(regressData)
```



If the data files are in Axona or Neuralynx format, they can be exported to HDF5 file

```
ndata.save_to_hdf5()
```

Datasets can be saved individually as well

```
spike.save_to_hdf5()
spat.save_to_hdf5()
```

Parametric results of all the analysis performed can be obtained by

```
results= ndata.get_results() # Returns the results in OrderedDict
print(results)
```

Results from individual data objects can also be retrieved similarly

```
spike_results = spike.get_results()
spat_results = spat.get_results()
```

## Analysis of head-directional cell

Change data filename/paths for the new unit similar to what was done for the place cell information Load new data and set the unit number. No need to reassign to ndata, as Python assignments are by reference, not by value.

```
ndata.set_data_format('NWB')
```

```
data_dir= \full\file\directory\of\head\directional\HDF5\data\'
spat.set_filename(data_dir + '120412_1.hdf5+/processing/Behavioural/Position')
spike.set_filename(data_dir + '120412_1.hdf5+/processing/Shank/3')
```

```
spat.load()
spike.load()
```

```
spike.set_unit_no(1)
```

Reset results to omit parametric output of previously analysed unit. This can be done before loading the new datasets or at any stage of the analysis.


```
ndata.reset_results()
```

Or, results can be reset using individual data objects

```
spat.reset_results()
spike.reset_results()
```

### Head-directional firing rate analysis and plot


```
hdData= ndata.hd_rate()
fig= nc_plot.hd_firing(hdData)
```



### Head directional shuffling analysis and plot


Number of bins for the histogram of the shuffled data is set to 100

```
hshuffleData= ndata.hd_shuffle(nshuff = 500, limit=0, bins= 100)
fig= nc_plot.hd_shuffle(hshuffleData)
```




## Head directional time-shift analysis

```
hshiftData= ndata.hd_shift(shift_ind=np.arange(-10, 10))
fig= nc_plot.hd_time_shift(hshiftData)
```



## Head directional multiple regression

```
regressData= ndata.multiple_regression()
fig= nc_plot.multiple_regression(regressData)
```



## Analysis of spike-train dynamics

### Changing the data filename/paths for the new unit

```
data_dir= \full\file\directory\of\HDF5\data\'
spat.set_filename(data_dir + '112512_1.hdf5+/processing/Behavioural/Position')
spike.set_filename(data_dir + '112512_1.hdf5+/processing/Shank/5')
```

```
spat.load()
spike.load()
```

```
spike.set_unit_no(1)
```


Reset results to omit parametric output of previously analysed unit

```
ndata.reset_results()
```

### Waveform properties of the unit

```
graphData= ndata.wave_property()
fig= nc_plot.wave_property(graphData, [int (spike.get_total_channels()/2), 2])
```






### Inter-spike interval (ISI) histogram

The number of bins for histogram is 350, and the maximum ISI to bin for is 350ms. This implies each bin represents 1msec interval. 'graphData' term will be used repeatedly from now on for reusing the memory


```
graphData= ndata.isi(bins = 350, bound = [0, 350])
fig= nc_plot.isi(graphData)
```



### ISI autocorrelation histogram for longer length

Binsize is 1msec, and autocorrelation is performed from -350ms to +350ms


```
graphData= ndata.isi_corr(bins = 1, bound = [-350, 350])
fig= nc_plot.isi_corr(graphData)
```



## ISI autocorrelation histogram for shorter length

Binsize is 1msec, and autocorrelation is performed from -10ms to +10ms


```
graphData= ndata.isi_corr(bins = 1, bound = [-10, 10])  
fig= nc_plot.isi_corr(graphData)
```



## Theta modulation Index analysis

Input parameters are for [Frequency,  $\tau_1$ ,  $\tau_2$ ] and provides the starting value, lower, and upper bound for the fitted sinusoidal equation. Binsize and temporal bound are that of ISI autocorrelation histogram

```
graphData= ndata.theta_index( start = [6, 0.1, 0.05], \  
                               lower = [4, 0, 0], \  
                               upper = [14, 5, 0.1], \  
                               bins = 1, bound = [-350, 350])  
fig= nc_plot.theta_cell(graphData)
```



Above analyses can also be done using the spike data itself as it does not require information from other data object. For example,

```
graphData= spike.isi(bins = 350, bound = [0, 350])
fig= nc_plot.isi(graphData)
```

## Analysis of rhythmicity of LFP and spike-to-LFP phase relationships

### Import NLfp class

```
from neurochat.nc_lfp import NLfp
```

**Instantiate LFP data object, set the filename, load data, and add to ndata**

```
lfp= NLfp(system= 'NWB')
```

```
lfp.set_filename(data_dir+ '\\112512_1.hdf5+/processing/Neural Continuous/LFP/eeg')
```

```
lfp.load()
```


```
ndata.lfp= lfp
```

### LFP frequency spectrum analysis

Hanning window of 2sec with 1sec overlap and number of FFT components= 2048. ptype is 'psd' which means power-spectral density. Other option can be 'power'. prefilter set 'True' for pre-filtering the LFP signal with a bandpass filter as set by filtset. filtset= [filter order, lower


cutoff frequency, higher cutoff frequency, type of filtering]. fmax defines the maximum frequency to analyse. db set to 'True' will convert the spectrogram in decibel unit. tr set to 'True' creates a time-resolved spectrogram with 'window'-resolution and 'overlap' amount of signal overlap. tr set to 'False' calculates the spectrogram using Welch's method. This function can also be similarly called as `ndata.spectrum()`

```
graphData= lfp.spectrum(window = 2, noverlap = 1, nfft = 2048, ptype = '
psd', \
    prefilt = True, filtset = [10, 1.5, 40, 'bandpass'], \
    fmax = 40, db = False, tr = False)
fig= nc_plot.lfp_spectrum(graphData)
```



After setting tr as True and db = True


```
graphData= lfp.spectrum(window = 2, noverlap = 1, nfft = 2048, ptype = '
psd', \
    prefilt = True, filtset = [10, 1.5, 40, 'bandpass'], \
    fmax = 40, db = True, tr = True)
fig= nc_plot.lfp_spectrum_tr(graphData)
```



## Spike-LFP phase distribution

$fwin = [6, 12]$  means that the phase of the spike are sought in the LFP band of 6Hz to 12 Hz. The minimum power of this band to be accepted to carry significant theta is 0.2 times the total LFP power, and that of the amplitude of the band signal is 0.15 times the amplitude of the LFP signal. The LFP signal is prefiltered using the filtset parameters.

```
graphData= ndata.phase_dist(binsize = 5, rbinsize = 2, fwin = [6, 12],\
                             pratio = 0.1, aratio = 0.15, filtset = [10, 1.5, 40, 'bandpass']
                             )
fig= nc_plot.spike_phase(graphData)
```



**The analysis can be performed from both the NLfp() and NSpike() objects**

Using the lfp object:

```
graphData= lfp.phase_dist(spike.get_unit_stamp(), binsize = 5, rbinsize  
= 2, fwin = [6, 12],\  
    pratio = 0.1, aratio = 0.15, filtset = [10, 1.5, 40, 'bandpass']  
)  
fig= nc_plot.spike_phase(graphData)
```


Using the spike object:

```
graphData= spike.phase_dist(lfp = lfp , binsize = 5, rbinsize = 2, fwin  
= [6, 12],\  
    pratio = 0.1, aratio = 0.15, filtset = [10, 1.5, 40, 'bandpass']  
)  
fig= nc_plot.spike_phase(graphData)
```

**Analysis of phase-locking value (PLV), spike-field coherence (SFC), and spike-triggered average (STA)**

Window of the LFP chunks in reference to the spike timestamps is set to -400ms to +400ms Frequency of interest for the analysis is set as 2Hz to 30Hz


```
graphData= ndata.plv(window = [-0.4, 0.4], fwin = [2, 40])  
fig= nc_plot.plv(graphData)
```



*Time-resolved* as set by mode= 'tr'. nsample implies number of randomly selected spikes around which the LFP signals are cut for phase-locking analysis slide gives the time in ms by which the window is shifted from left to right to obtain the time-resolved phase-locking analysis

```
graphData= ndata.plv(window = [-0.4, 0.4], nfft = 1024, mode = 'tr', nsample = 2000, slide = 25, fwin = [2, 40])
fig= nc_plot.plv_tr(graphData)
```





In most of the cases where composite information are required and `ndata` is not used, the spike timestamp is provided as the first argument to the methods followed by other information. Because, in such cases only information required by the analysis from the spike object is the timestamps of individual spikes in the train. For example,

```
graphData= ndata.plv(window = [-0.4, 0.4], fwin = [2, 40])
fig= nc_plot.plv(graphData)
```

gives the same result as the codes given below:

```
graphData= lfp.plv(spike.get_unit_stamp(), window = [-0.4, 0.4], fwin =
[2, 30])
fig= nc_plot.plv(graphData)
```

## Use of Nhdf class

### Import and instantiate Nhdf class

```
from neurochat.nc_hdf import Nhdf
hdf= Nhdf()
```

### Store data using Nhdf object

Nhdf() resolves the filename and the path for storage of the data using

Nhdf().resolve\_pathname(data=data\_obj) where data\_obj can be a NSpatial(), NSpike() or NLfp() object

```
hdf.save_spatial(spat)
hdf.save_spike(spike)
hdf.save_lfp(lfp)
```

### This can also be done using

```
hdf.save_object(obj = spat)
hdf.save_object(obj = spike)
hdf.save_object(obj = lfp)
```

## Graphical data from individual analysis can be stored using the following codes

*path* is the path inside HDF5 file. Analysis data are always recommended to store in the */analysis/* path. But analysis for each unit+lfp pair is stored in one path under which graphical data from individual analyses are stored. The unique unit ID is established using the name resolving method Nhdf().resolve\_analysis\_path() which utilizes the filename of the recorded data, electrode/tetrode number, eeg channel ID and the unit number. *name* is the name of the analysis following the unit ID, i.e. 'plv' etc. *graph\_data* are the dictionary data that are plotted using the functions in nc\_plot

```

unit_id= hdf.resolve_analysis_path(spike = spike, lfp = lfp)

hdf_name= hdf.resolve_hdfname(data=spike) # Resolve HDF5 filename
hdf.set_filename(hdf_name) # NeuroChaT opens the file as file-object as soon as new filename is set.

print(unit_id)

hdf.save_dict_recursive(path = '/analysis/' + unit_id+ '/',
                        name = 'plv', data = graphData)

```

### Analysis results can be stored by

```

results= ndata.get_results()
hdf.save_dict_recursive(path = '/analysis/' + unit_id+ '/',
                        name = 'results', data = results)

```

### Apart from that data and attributes to any group or dataset can be added using

Set `create_group` to 'True' it will create the path if does not already exist

```

hdf.save_dataset(path = '/path/to/group/', name = 'name_of_dataset', data = data_to_store, create_group = True)
hdf.save_attributes(path= '/path/to/group/or/dataset/', attr = dict_of_attributes)

```

## Use of NeuroChaT class

### Import NeuroChaT class and instantiate

```

from neurochat.nc_control import NeuroChaT
nc= NeuroChaT()

```

### Convert files in Axona format to NWB files specified in an Excel list

```

excel_file= '\\full\\file\\name\\of\\Excel\\list.xlsx'
nc.convert_to_nwb(excel_file)

```

### Verify units provided in an Excel list before batch-mode analysis

```

excel_file= '\\full\\file\\name\\of\\Excel\\list.xlsx.xlsx'
nc.verify_units(excel_file)

```

## Evaluate the quality of clustering from a list provided in an Excel file

```
excel_file= '\full\file\name\of\Excel\list.xlsx.xlsx'  
nc.cluster_evaluate(excel_file)
```

## Evaluate similarity of clusters

The excel list contains paired list of units to be compared for similarity

```
excel_file= 'C:\\Users\\Raju\\Google Drive\\Sample Data for NC\\Comparis  
on results_from NeuroChaT_pawels_data.xlsx'  
nc.cluster_evaluate(excel_file)
```

## Analysis using NeuroChaT

Analysis using NeuroChaT class is always done with the help of Configuration class where the user specifies all the data, intended analyses, input parameters etc.

### Configuration class

Import, instantiate, set the filename and load configuration from the file. This class uses nc\_defaults.py module for importing default analyses and parameters.

```
from neurochat.nc_config import Configuration  
  
config= Configuration()  
  
config.set_config_file('\full\file\name\of\grid_config.ncfg')  
  
config.load_config()
```

Set configuration to NeuroChaT object

```
nc.set_configuration(config)
```

Start analysis. This will 'read' the instructions from the config object and execute accordingly

```
nc.start()
```

Use *get\_* and *set\_* functions also known as getters and setters for accessing and setting values of interest. For example, \* Getting and setting parameters:

```

param_list= config.get_param_list() # List of all parameters as dictionary keys

params_by_analysis= config.get_params_by_analysis(analysis= 'isi')
print(params_by_analysis)

param_val= config.get_params(name = 'isi_length') # name is the list of parameters or the name of a single parameter'
print(param_val)

config.set_param(name = 'isi_bin', value = 2)

```

### Getting and setting analyses

```

list_of_analyses= config.get_analysis_list() # List of all analysis
print(list_of_analyses)

analysis_checked= config.get_analysis(name = 'isi') # If 'True', analysis is set to be done
print(analysis_checked)

config.set_analysis(name = 'theta_skip_cell', value = False) # Analysis of theta skipping cell turned off

```

Analyses can be performed in different modes, namely: 1. 'Single Unit'-one cell at time, value '0' 2. 'Single Session'- all the cells in one recording at a time, value '1' 3. 'Listed Units'- all the cells listed in one Excel file, value '2'

### Getting and setting analysis mode

```

print(config.get_analysis_mode())

config.set_analysis_mode(analysis_mode = 'Single Unit') # Can also set analysis_mode = 0

```

What type of data file need to be specified depends on the type of mode and the format of the data Please refer to the Configuration class for more such methods. Here, we show an example of settingh Axona data and an example of batch mode analysis

### Specifying Axona files for analyses

```

data_dir= '\path\to\recorded\Axona\data\'

config.set_analysis_mode(0) # For 'Single Unit' analysis

config.set_spatial_file(spatial_file = data_dir+ '040513_1_1.txt')
config.set_spike_file(spike_file = data_dir + '040513_1.6')

config.set_unit_no(3)

```

We are interested in only certain analyses. So, we first turn off all the analyses:

```
config.set_analysis(name = 'all', value = False) # 'all' for setting all the analyses
```

### Specify new analyses

```
config.set_analysis(name = ['loc_rate', 'loc_shuffle', 'loc_time_lapse'], value = True) # See nc_defaults for names of the analyses
```

Let us use default parameters for ease of understanding. NeuroChaT() always saves the graphics in a file. Let us set the file in 'PDF' or 'pdf' format. Other option is 'Postscript' or 'ps'

```
config.set_graphic_format(graphic_format = 'PDF')
```

### Set this configuration for NeuroChaT's use

```
nc.set_configuration(config)
```

Save this configuration to a file for future use. This file can be edited using any standard text-editing software

```
config.save_config('\full\file\name\of\place_config.ncfg')
```

Once the configuration file is set to NeuroChaT object, all of its methods can be used by NeuroChaT itself. For example, the configuration can be loaded from and saved to file using the NeuroChaT object. It works this way- if NeuroChaT cannot find a method within itself, it at first searches in the Configuration object. If not found, it looks into composing object NData() for the function. This process is called delegation. The precedence for delegation is Configuration() > NData()

```
nc.set_config_file('\full\file\name\of\place_config.ncfg')
nc.load_config()
```

```
nc.set_analysis_mode(0) # Analysis mode set to 'Single Unit' in Configuration object
```

Once the analyses are done, NeuroChaT saves the pdf in respective data folder. It always stores the NWB-converted file if the latter does not exist and stores the graphics data and the parametric results in the files. Along with that, parametric results and names of output PDF and NWB files can be obtained by using following codes which return them in Pandas DataFrame.

```

results_df= nc.get_results()
print(results_df)
output_filename_df= nc.get_output_files()
print(output_filename_df)

TT6_SS_4_eeg 9.730225 23.065766 21.465309 241.153945
Mean amplitude Std width Mean height Theta Index
TT6_SS_4_eeg 203.199722 64.520798 204.495651 0.714889
TI fit freq Hz TI fit tau1 sec ... Mult Rsq
TT6_SS_4_eeg 8.808084 0.229588 ... 0.222366
Semi Rsq Loc Semi Rsq HD Semi Rsq Speed Semi Rsq Ang Vel
TT6_SS_4_eeg 0.15583 0.002322 0.03563 0.001138
Semi Rsq Dist Border DR HP DR SP DR AP DR BP
TT6_SS_4_eeg 0.001403 0.085843 0.340246 0.190116 0.159364
[1 rows x 88 columns]
Graphics Files
TT6_SS_4_eeg C:\Users\Raju\Google Drive\Sample Data for NC...
NWB Files
TT6_SS_4_eeg C:\Users\Raju\Google Drive\Sample Data for NC...

```

These files can be exported for future use using DataFrame's io utilities:

```

import pandas as pd
writer= pd.ExcelWriter('\full\file\path\to\parametric_results.xlsx') # s
et-up writing engine
results_df.to_excel(writer, 'Sheet1') # write to file
output_filename_df.to_excel(writer, 'Sheet2')

```

While the graphical interface provides an easier means for performing almost all of the abovementioned functionalities, NeuroChaT and its constituent classes works as the 'engine' behind those tasks.

## Use NClust class

### Import and instantiate NClust

Although we are initializing it with already defined spike object, we could similarly set the filename and unit and load the composing spike object as we do for any other spike object itself NClust also performs some of the analysis that spike object does, i.e. analysing waveform properties, ISI histogram, PSTH etc. See nc\_clust.py module to learn more about this aspect.

```
from neurochat.nc_clust import NClust
clust= NClust(spike= spike)
```

This object is intended for facilitating analysis pertaining to clustering algorithm and cluster quality measurements. Following are some of the example methods:

#### **Remove null channels if any**

```
off_chan= clust.remove_null_chan()
```

#### **Resample wave by intended factor**

```
wave, time= clust.resample_wave(factor= 2) # Resampling factor is 2
```

#### **Align waves by peaks for better estimation of waveform features**

```
clust.align_wave_peak()
aligned_wave= clust.getWaveform()
```

#### **Get the channel with highest waveform energy, peak at the channel , and the index of the peak**

```
peak, peak_chan, maxInd= clust.get_max_wave_chan()
```

#### **Get the Principle Components of the waveforms**

```
pc= clust.get_wave_pc(npc = 2) # 2 PC in each channel
print(pc)
```

#### **Get features for clustering**

```
feat= clust.get_feat(npc = 2) # Consist of waveform peaks, troughs and 2 PC components in each channel
```

#### **Get fetures of clustered units**

```
unit_feat= clust.get_feat_by_unit(unit_no = 3)
```

#### **Get waveforms by unit number**

```
waves= clust.get_unit_waves()
```

#### **Clustering quality evaluation**



If unit\_no set to 0 all units are evaluated with a matrix output for pairwise comparison. Otherwise, maximum Bhattacharyya distance (BC) and minimum Hellinger distance (Dh) for the specified unit are returned

```
bc, dh = clust.cluster_separation(unit_no = 0)
```

### **Evaluationg unit similarity**

```
clust_1 = NClust()
clust_1.load(filename = '\\full\\file\\directory\\of\\spike\\data_1', system =
'NWB') # An alternative approach for loading spike data

clust_2 = NClust()
clust_2.load(filename = '\\fullfile\\directory\\of\\spike\\data_2', system =
'NWB') # An alternative approach for loading spike data

bc, dh = clust_1.cluster_similarity(nclust= clust_2, unit_1= 3, unit_2=
3) # unit_1 and unit_2 are the comparable units
```

## C. Input parameter description

### Waveform properties

NO PARAMETER

### Inter-spike interval (ISI)

isiBin	2	1-100	ms	Bin size of the ISI histogram
isiLength	350	10-1000	ms	Length of ISI histogram
isiLogBins	70	10-100		
isiLogLength	350	10-1000	ms	

### ISI Autocorrelation

isiCorrBinSh	1	1-10	ms	Bin size of the ISI correlation histogram obtained on short lags
isiCorrLenSh	10	5-50	ms	Length of the ISI correlation histogram obtained on short lags
isiCorrBinLong	2	1-50	ms	Bin size of the ISI correlation histogram obtained on long lags
isiCorrLenL	350	10-1000	ms	Length of the ISI correlation histogram obtained on long lags

### Theta-modulated Cell Index

thetaCellFreqMin	6	1-10	Hz	Lower limit of the theta band frequency in curve fitting
------------------	---	------	----	--

thetaCellFreqMax	12	8-16	Hz	Upper limit of the theta band frequency in curve fitting
thetaCellFreqStart	6	5-10	Hz	Starting value of the theta band frequency in curve fitting
thetaCellTaulMax	5	0.5-10	sec	Upper limit of the decay constant $\tau_1$ decay constant in curve fitting
thetaCellTaulStart	0.1	0-15	sec	Starting value of the decay constant $\tau_1$ decay constant in curve fitting
thetaCellTau2Max	0.05	0-0.1	sec	Upper limit of the decay constant $\tau_2$ decay constant in curve fitting
thetaCellTau2Start	0.05	0-0.1	sec	Starting value of the decay constant $\tau_2$ decay constant in curve fitting

### Theta-skipping Cell Index

NO PARAMETER	Parameters from theta-modulated cell index are be used
--------------	--

### Burst Property

burstThresh	5	1-15	ms	Minimum ISI between consecutive spikes in a burst
spikesToBurst	2	2-10	ms	Minimum number of consecutive spikes with burstThresh for a burst
ibiThresh	50	5-1000	ms	Minimum inter-burst interval between two bursting groups of spikes

### Spike Rate vs Running Speed

speedBin	1	1-10	cm/sec	Size of the speed bin for histogram
speedMin	0	0-10	cm/sec	Minimum acceptable speed to analyse
speedMax	40	10-200	cm/sec	Maximum limit on the speed to analyse
speedKernLen	3	1-25, odds	samples	Length of moving-average smoothing kernel of recorded speed
speedRateKernLen	3	1-7, odds	bins	Length of moving-average smoothing kernel of the spike rate

### Spike Rate vs Angular Velocity

angVelBin	10	1-50	deg/sec	Size of the angular velocity bin for histogram
angVelMin	-200	-500-0	deg/sec	Minimum acceptable angular velocity to analyse
angVelMax	200	0-500	deg/sec	Maximum limit on the angular velocity to analyse
angVelCutoff	10	0-100	deg/sec	
angVelKernLen	3	1-25, odds	samples	Length of moving-average smoothing kernel of calculated angular velocity
angVelRateKernLen	3	1-5, odds	bins	Length of moving-average smoothing kernel of the spike rate

## Spike Rate vs Head Direction

hdBin	5	factors of 360	degree	Size of the head directional bin for histogram
hdAngVelCutoff	30	0:5:100	deg/sec	Lower limit of the acceptable angular velocity for avoiding noise from jerking of the head
hdRateKernLen	5	1-11	bins	Length of moving-average smoothing kernel of the spike rate

## Head Directional Shuffling Analysis

hdShuffleTotal	500	100-10000		Number of shuffles for head directional shuffling analysis
hdShuffleLimit	0	0:2:500	sec	Upper limit of the shuffled shifted time of spikes
hdShuffleNoBins	100	10:10:200	bins	Number of bins for displaying the distribution of specificity measures

## Head Directional Time Lapse Analysis

NO PARAMETER				
--------------	--	--	--	--

## Head Directional Time Shift Analysis

hdShiftMax	10	1-100	indices	Maximum number of spatial samples by which spike-events are shifted forward
------------	----	-------	---------	---

hdShiftMin	-10	-1 to -100	indices	Maximum number of spatial samples by which spike-events are shifted backward
hdShiftStep	1	1-3		Steps of samples spikes are shifted with

### Spike Rate vs Location

locPixelSize	3	1-100	cm	Pixel size (bin size) for spatial firing 2D histogram
locChopBound	5	3-20	pixels	Upper limit on number of empty rows and columns which are chopped off from the firing rate map
locRateFilter	Box	Box/Gaussian		Type of smoothing kernel of the firing rate map
locRateKernLen	5	1-11, odds	pixels/ no unit	Number of pixels for box filter and standard deviation for Gaussian filter

### Locational Shuffling Analysis

locShuffleTotal	500	100-10000		Number of shuffles for locational shuffling analysis
locShuffleLimit	0	0:2:500	sec	Upper limit of the shuffled shifted time of spikes
locShuffleNoBins	100	10:10:200	bins	Number of bins for displaying the distribution of specificity measures

## Locational Time Lapse Analysis

NO PARAMETER

## Locational Time Shift Analysis

locShiftMax	10	1-100	indices	Maximum number of spatial samples by which spike-events are shifted forward
locShiftMin	-10	-1 to -100	indices	Maximum number of spatial samples by which spike-events are shifted backward
locShiftStep	1	1-3		Steps of samples spikes are shifted with

## Spatial Autocorrelation

spatialCorrMinObs	20	1-100	pixels	Minumum number of overlapping pixels between original and the shifted firing rate map
rotCorrBin	3	factors of 360	degree	Steps of firing rate map rotation
spatialCorrFilter	Box	Box/Gaussian		Type of smoothing kernel of the firing rate map
spatialCorrKernLen	5	1-11, odds		Number of pixels for box filter and standard deviation for Gaussian filter

## Grid Cell Analysis

Parameters from spatial autocorrelation are also used				
gridAngTol	2	1-5	degree	Accepted differences between angular measures to consider them same
gridAngBin	3	Factors of 360 less than 45	degree	Size of angular bins to measure rotational correlation

## Border Cell Analysis

borderFiringThresh	0.1	0:0.05:1		Minimum firing rate to maximum firing rate ratio to define as active pixels
borderAngBin	3	Factors of 360 less than 45		Size of angular bin for circular-linear firing rate map
borderStairSteps	5	4-10		Number for steps for stair plot of border firing

## Gradient Cell Analysis

gradAsympLim	0.25	0.1:0.05:1		Range of asymptotic parameter 'a' on Goempertz function, e.g., $\pm 0.25$
gradDisplaceLim	0.25	0.1:0.05:1		Range of displacement parameter 'b' on Goempertz function, e.g., $\pm 0.25$
gradGrowthRateLim	0.5	0.1:0.05:1		Range of growth rate parameter 'c' on Goempertz function, e.g., $\pm 0.5$



## Multiple Regression

multiRegInterval	0.1	0.1:0.1:1	sec	Interval of spatial samples for multiple regression
multiRegEpisode	120	60:30:300	sec	Duration for each replication of multiple regression
multiRegNoRep	1000	100:100:2000		Number of replications for multiple regression

## Interdependence Analysis

NO PARAMETER	Parameters from other analysis should be used
--------------	---

## LFP Frequency Spectrum

lfpPreFiltLowCut	1.5	0.1:0.1:4	Hz	LFP prefiltering lower cutoff frequency
lfpPreFiltHighCut	40	10:5:500	Hz	LFP prefiltering higher cutoff frequency
lfpPreFiltOrder	5	1-20		LFP prefiltering butterworth filter order
lfpPwelchSegSize	2	0.5:0.5:100	sec	LFP segment size for Welch's method for spectral density estimation
lfpPwelchOverlap	1	0.5:0.5:50	sec	Overlap between LFP segments for Welch's method for spectral density estimation
lfpPwelchNfft	1024	128:128:8192		NFFT for Welch's method for spectral density estimation
lfpPwelchFreqMax	40	10:5:500	Hz	Maximum frequency to display

lfpStftSegSize	2	0.5:0.5:100	sec	LFP segment size for short-time Fourier transform
lfpStftOverlap	1	0.5:0.5:50	sec	Overlap between LFP segments for short-time Fourier transform
lfpStftNfft	1024	128:128:8192		NFFT for short-time Fourier transform
lfpStftFreqMax	40	10:5:500	Hz	Maximum frequency to display

### Unit LFP-Phase Distribution

phaseFreqMin	6	1-10	Hz	Lower frequency of LFP band for analysis of phase locking
phaseFreqMax	12	1-10	Hz	Higher frequency of LFP band for analysis of phase locking
phasePowerThresh	0.1	0:0.05:1		Minimum band power to overall power of acceptable LFP segments
phaseAmpThresh	0.15	0:0.05:1		Minimum segment amplitude to overall amplitude of acceptable LFP segments
phaseBin	5	Factors of 360	degree	Size of phase bins for circular histogram of spike-phases
phaseRasterBin	2	1-15	degree	Size of phase bins for raster of spike-phases

### Unit LFP-Phase Locking

phaseLockWinLow	-0.4	-1:0.05:-0.1	sec	Lower limit of the LFP segments in reference to spike-events
phaseLockWinUp	0.4	0.1:0.05:1	sec	Upper limit of the LFP segments in reference to spike-events
phaseLockNfft	1024	128:128:8192		NFFT for Fourier transform

phaseLockFreqMax	40	10:5:500	Hz	Maximum frequency to analyse and display
------------------	----	----------	----	--

## **D. Output parameter description**

### **Waveform Properties**

Mean Spiking Freq= Spiking frequency over the entire trial

Mean amplitude= mean of amplitude. Amplitude is measured by the difference between the first negative peak and first positive peak,

Std amplitude= Standard deviation of amplitude

Mean height= mean of height. Height is measured by the difference between the first positive peak and the overall minimum value of the spike

Std height= Standard deviation of height

Mean width= Mean of spike width taken at 25% of the amplitude

Std width= Standard deviation of the width

### **Inter-spike Interval (ISI)**

No output parameter

### **ISI Autocorrelation**

No output parameter

### **Theta-modulated Cell Index**

Theta Index= Ratio of the sinusoid amplitude and the constant term in the model

TI fit freq Hz= Theta frequency from the model fitting in Hz

TI fit taul sec= Decay constant for the exponential modulation of the sinusoidal component in the model

TI adj Rsq= Goodness of fit of the model for the ISI autocorrelation data

TI Pearse R= Goodness of fit Pearson's R between the original and model-fit values for the ISI autocorrelation data

TI Pearse P= Goodness of fit Pearson's P between the original and model-fit values for the ISI autocorrelation data

### **Theta-skipping Cell Index**

Theta Skip Index= Theta Skipping Index

TS jump factor= Ratio of the amplitude of the theta and delta band sinusoids

TS f1 freq Hz= Fitted frequency of the faster sinusoid in Hz

TS f2 freq Hz= Fitted frequency of the slower sinusoid,

TS freq ratio= Ratio of f1 to f2

TS taul Hz= Decay constant for the exponential modulation of the sinusoidal component in the model

TS adj Rsq= Goodness of fit of the model for the ISI autocorrelation data

TS Pearse R= Goodness of fit Pearson's R between the original and model-fit values for the ISI autocorrelation data

TS Pearse P= Goodness of fit Pearson's P between the original and model-fit values for the ISI autocorrelation data

### **Burst Property**

Total burst= Total number of bursts calculated

Total bursting spikes= Total number of spikes constituting bursts

Mean bursting ISI ms= Mean inter-spike-interval for the bursting spikes only (ms)

Std bursting ISI ms = Standard deviation of inter-spike-interval for the bursting spikes only (ms)

Mean spikes per burst= Average number of spikes in bursts

Std spikes per burst= Standard deviation of the number of spikes in bursts

Mean burst duration= Mean of the duration of bursts in ms

Std burst duration= Standard deviation of the duration of bursts in ms

Mean duty cycle= Mean of the duty cycles. Duty cycle is the portion of the inter-burst interval during which the burst fires (burst duration/ inter-burst interval)

Std duty cycle= Standard deviation of the duty cycles

Mean IBI= Mean of inter-burst Intervals in ms

Std IBI= Standard deviation of inter-burst intervals in ms

Propensity to burst= Total bursting spikes/total spikes in the cluster

### **Spike Rate vs Running Speed**

Speed Skaggs= Skaggs information content for speed vs spiking events in bits/sec

Speed Pears R= Goodness of fit Pearson's R between spike firing rate at different speeds and fitted straight line

Speed Pears P= Goodness of fit Pearson's P between spike firing rate at different speeds and fitted straight line

## **Spike Rate vs Angular Velocity**

Ang Vel Left Pears R= Goodness of fit Pearson's R between spike firing rate at different counter-clockwise (-ve) angular head velocity and fitted straight line

Ang Vel Left Pears P= Goodness of fit Pearson's P between spike firing rate at different counter-clockwise (-ve) angular head velocity and fitted straight line

Ang Vel Right Pears R= Goodness of fit Pearson's R between spike firing rate at different clockwise (+ve) angular head velocity and fitted straight line

Ang Vel Right Pears P= Goodness of fit Pearson's P between spike firing rate at different clockwise (+ve) angular head velocity and fitted straight line

## **Spike Rate vs Head Direction**

HD Skaggs= Head directional Skaggs information content

HD Rayl Z= Rayleigh Z for the head-directional firing rate

HD Rayl P= Rayleigh P for the head directional firing rate

HD von Mises K= von Mises concentration parameter  $\kappa$

HD Mean= Vector mean or preferred head direction (degree) of the unit to fire

HD Mean Rate= Firing rate in preferred direction

HD Res Vect= Resultant vector length of head-directional firing rate

HD Peak Rate= Peak firing rate in the head-directional tuning curve

HD Peak= Head-direction at which peak firing rate occurs

HD Half Width= Width of the tuning curve measured at 50% of the peak firing rate(degree)

HD Peak CW= Peak firing direction during clockwise head-directional movement

HD Peak CCW= Peak firing direction during counter-clockwise head-directional movement

HD Peak Rate CW= Peak firing rate during clockwise head-directional movement

HD Peak Rate CCW= Peak firing direction during counter-clockwise head-directional movement

HD Delta= Separation angle between peak firing direction during clockwise and counter-clockwise head-directional movement

### **Head Directional Shuffling Analysis**

HD Shuff Rayl Z Per 95= 95<sup>th</sup> percentile of the distribution of Rayleigh Z parameter for the head-directional tuning curves obtained from shuffling of spike-events

HD Shuff von Mises K Per 95= 95<sup>th</sup> percentile of the distribution of von Mises concentration parameter  $\kappa$  for the head-directional tuning curves obtained from shuffling of spike-events

### **Head Directional Time Lapse Analysis**

No Parameter

### **Head Directional Time Shift Analysis**

HD ATI= Anticipatory time interval for the head-directional cells, measured as the time-shift where the counter-/clockwise head-directions are same, or the separation angle becomes zero.

HD Opt Shift Skaggs= Time shift which maximizes the information content in HD tuning of spiking events

HD Opt Shift Peak Rate= Time shift which maximizes the peak firing rate in HD tuning of spiking events

### **Spike Rate vs Location**

Spatial Skaggs= Information content of spatial firing map

Spatial Sparsity= The fraction of the environment in which the cell is active (max=1, min =0)



Spatial Coherence= Measure of orderliness of the local firing pattern (max=1. Min = -1). Or, simply the correlation between raw firing map, and smoothed firing map (value at each pixel is replaced the 8 neighbouring pixels of the non-smooth map). See Muller & Kubie 1989

### **Locational Shuffling Analysis**

Loc Skaggs 95= 95<sup>th</sup> percentile of the distribution of Skaggs information content for the firing rate map obtained from shuffling of spike-events

Loc Sparsity 05= 95<sup>th</sup> percentile of the distribution of sparsity for the firing rate map obtained from shuffling of spike-events

Loc Coherence 95= 95<sup>th</sup> percentile of the distribution of coherence for the firing rate map obtained from shuffling of spike-events

### **Locational Time Shift Analysis**

Loc Opt Shift Skaggs= Time shift which maximizes the Skaggs information content in spatial firing map

Loc Opt Shift Sparsity= Time shift which minimizes the sparsity in spatial firing map

Loc Opt Shift Coherence= Time shift which maximizes the spatial coherence in spatial firing map

### **Spatial Autocorrelation**

No Parameter

### **Grid Cell Analysis**

Is Grid= Indicates if the unit is a Grid cell or not (1= yes, 0= no)

Grid Mean Alpha= Average of the angles each arm of the hexagon (formed from the peaks of the firing fields) forms with the centre of the spatial autocorrelation

Grid Mean Psi= Mean angle between the arms of the central hexagon (formed from the peaks of the firing fields) in spatial autocorrelation

Grid Spacing= Average spacing between the peak firing fields forming the grid; Obtained from spatial autocorrelation

Grid Score= Gridness score as in DOI: 10.1126/science.1125572

Grid Orientation= Inclination of the central hexagonal firing field patterns with the X-axis

### **Border Cell Analysis**

Border Skaggs= Skaggs information content for border vs spike-rate

Border Ang Ext= Largest angular segment with non-zero histogram count in active pixel (>20% of maximum firing rate) vs angular distance histogram

### **Gradient cell Analysis**

Grad Pearse R= Goodness of fit Pearson's R between the calculated and model-fit rate of firing vs distance from border

Grad Pearse P= Goodness of fit Pearson's P between the calculated and model-fit rate of firing vs distance from border

Grad adj Rsq= Goodness of adjusted  $R^2$  between the calculated and model-fit rate of firing vs distance from border

Grad Max Growth Rate= Maximum rate of growth in firing rate in the fitted Goempertz function

Grad Inflect Dist= Distance from border where the growth of firing rate is maximum

## **Multiple Regression**

Multi Rsq= Goodness of fit of the linear equation with the observed spike rate. Alternatively, it is a measure of the amount of variance explained in firing rate by all the independent variables

Semi Rsq Loc= Explained variance in firing rate by the location alone

Semi Rsq HD= Explained variance in firing rate by the head direction alone

Semi Rsq Speed= Explained variance in firing rate by the running speed alone

Semi Rsq Ang Vel= Explained variance in firing rate by the angular velocity alone

Semi Rsq Dist Border= Explained variance in firing rate by the border as a variable alone

## **Interdependence Analysis**

DR HP= Distributive ratio for predicting head direction tuning curve (H) from spatial firing map (P)

DR SP= Distributive ratio for predicting spike rate vs running speed curve (S) from spatial firing map (P)

DR AP= Distributive ratio for predicting spike rate vs Angular velocity curve (A) from spatial firing map (P)

DR BP= Distributive ratio for predicting spike rate vs distance from border curve (B) from spatial firing map (P)

## **LFP Frequency Spectrum**

No parameter

## **Unit LFP-phase Distribution**

LFP Spike Mean Phase= Average LFP phase of the spikes

LFP Spike Mean Phase= Average no of spikes with Mean Phase

LFP Spike Phase Res Vect= Resultant vector on the distribution of spike-phases on LFP waves

## **Unit LFP-phase Locking**

No parameter